



OMI PADS PEMF Beyond mattress



Advanced Functionality Full body PEMF mat
for the overall well-being of the body

Description

High quality, easily foldable mattress with state-of-the-art, built-in battery-operated wireless controller unit. The mattress can be operated via the OMI Beyond app on any iOS, or Android device. The apps are free to download from the iOS Appstore, or the Android Play Store. It features over 100 preset programs. It works in the frequency range from 1 to 99 Hz with the most effective ranges from 3 to 40 Hz. Mattress size: 165 x 65 cm. Increased magnetic field intensity up to 1,000+ μ T. Adjustable Rectangular, Sine and Chainsaw Waveforms. Powered by DC power supply.

Recommended for

The mattress is recommended for people who live stressful life, feel less energetic and wish to improve their overall well-being.

Claims

The pulsing magnetic field promotes the body's own regeneration processes and helps

- ✓ to stimulate micro-circulation and cell metabolism
- ✓ to maximize mental clarity, focus and productivity
- ✓ to increase energy and vitality
- ✓ to aid relaxation
- ✓ to sleep better and have better general mood

Advantages

State-of-the-Art technology, advanced functionality, ergonomic and foldaway design, for both professional and home use. Completely safe, not painful to use, and there are no known side effects.

Instructions for use

Easy and convenient to use. No special training is required. Just download the OMI Beyond controller app from the iOS Appstore, or the android Play Store. Connect the controller unit to the mat and switch on the controller! All is very straightforward.

Open the App on your phone and follow the guides! You may use the OMI Beyond mat right out of the box.

For detailed use instructions, please visit the following page:

https://www.oxfordmedicals.com/How-to-use-the-OMI-PEMF-Beyond-Mat_b_53.html

How it works

Pulsed electromagnetic fields produced by the device induce electrical current inside the body that helps to restore the blocked ion flow and stimulates the body's own regeneration processes.

As a result of using the mattress the body attains a higher energy level, the fatigue is alleviated, helping to achieve improved overall well-being of the whole body.

Packaging contains: Mattress, operating unit, wall adapter, user manual

Documentation and others

Declaration of Conformity (LVD, EMC)

Wholesale Inquiries

All level of wholesale inquiries are welcome. To discuss further details, please [contact us](#).

Payment Terms: on delivery

Pick-up Point: Budapest

Minimum Order Quantity: [ask for quote](#)

Country of Origin: European Community

Delivery Term, Currency and Price: [ask for quote](#)

Package Details

Artificial leather carry-on bag in a mail order carton



Product Dimensions (cm): 81 x 43 x 9 cm

Product Net Weight (kg): 2.3 kg

Product Gross Weight (kg): 3.5 kg

Master Carton Dimensions (cm): 83 x 45 x 47 cm

Master Carton Gross Weight (kg): 18.5 kg

Master Carton Quantity (pc): 5 pcs

The product is not a medical device; it is not intended to diagnose, treat, cure or prevent any disease.

product category: wellness product

website: www.oxfordmedicals.com



Manufactured by
Oxford Medical Instruments Limited
2nd Floor Curzon House, 24 High Street, Banstead, Surrey SM7 2LJ, England

www.oxfordmedicals.com